


- 1) Add header on binary file with
 - 1) ID of the digitizer
 - 2) Data structure (in Setup directory?)
 - 3) Run number
- 2) Structure of the data director with separate subdirectory for tracker, PID, gamma e setup
- 3) The macro that convert binary to root file take the link between digitizer and anode section from separate configuration file. The macro should provide for different level of
 - min level all the original data are kept to have maximum flexibility during test.
 - max level only the useful data are kept and all the other are not copied to save space
- 4) Separation of the flag variable in two different variable (to be aligned with the digitizers documentation)
- 5)





DEFINITIONS

- *Hit / entry*: All the information relate to the signal in a single pad
entry is related to the organization of this information on a root or binry file,
hit is the corresponding physical phenomeno
- *Event*: It is the totality of hit generated by a single physical event (e.g.
charged particle crossing the detector)
- *Cluster*: Group of entries that are close in space and time (not necessarily
part of a physical event)
- *Row*: is the totality of pad or strip tha have same z-coordinate (in the final
detector a row is made or of 64x). It is divided in four segments correspdng to
the four section of the anode.
- *Strip*: Is the sensible area of the anode between two rows. Each strip is made of
four long segment (300 mm) each one corresponding to a section of the anode.

Inside the code *strip* and *row* are both called row. The proper raw have an index that run from 0 to 4 and 64x4 pads, the strips instead have an index between 5 and 10 and just 4 pads.


pad

strip segment



row

Strip (row inside code)



Full anode (group of all 4 anodes)

Single anode (single anodic board)

Anode section (group of pad and strip linked to a digitizer)
At each anode section is linked one and only one digitizer and a preamplifier

Anode section	row	anode board	ID dig	ID preampl
1	0	1	xxx	yyyy
2	1	1	xyy	yyyx
3	2	1
4	3	1
5	4	1
6	0	2
7	1	2
20	4	4		

Template for the anode map file.


Software nomenclature

Top view / External view of a single anodic board

The diagram shows a 600x60 grid representing a beam cross-section. The grid is divided into 6 horizontal rows, each 100 units high. The rows are labeled 0 to 5 on the left and 5 to 10 on the right. The total width is 600 units. The grid is composed of 60 columns, each 10 units wide. The first 59 columns are labeled '59' on the left. The last 5 columns of each row are labeled '4', '3', '2', '1', and '0' on the right. A blue arrow labeled 'beam' points upwards from the bottom center of the grid.

Electronics nomenclature

Bottom view / internal view of a single anodic board



Electronics nomenclature

Bottom view / internal view of a single anodic board

row

View from Bottom side (Scale 1:1.08154877784988)

E

2

D

C

B₃

A

PI6


PI5

PI4


PI3

PI2

PI1



beam



JA	row A + PI1	row 0 + row 5
JB	row B + PI2	row 1 + row 6
JC	row C + PI3	row 2 + row 7
JD	row D + PI4	row 3 + row 8
JE	row E + PI5 + PI6	row 4 + row 9 + row 10


Top view / External view of the full anode with an highlight (red) on a single section of the anode

[illegible]

I.N.F.N. Napoli
AnodoNumenV2022.PcbDoc
cod. Numen ver 1
04/04/2022

Bot Poste

DRILL DRAWING
VIEW TOP




The active area of the detector
start from +2 (To be verified)

Pinout of different elements

Pad number	output channel PA
0	1
1	2
2	3
3	4
⋮	⋮
59	60
Strip	61
Strip2	62 (if row is 4)


Pin number SAMTEC ERF8	channel PA
1	bias 0:31
2	bias 0:31
3	N.C
4	N.C
5	16
6	0
7	17
8	1
9	18
10	2
⋮	⋮
33	30
34	14
35	31
36	15
37	N.C.
38	N.C.
39	GND
40	GND
41	GND
42	GND
43	N.C.
44	N.C.
45	48
46	32
47	49
48	33
49	50
50	34
⋮	⋮
71	61
72	45
73	62
74	46
75	63
76	47
77	N.C.
78	N.C.
79	bias 32:63
80	bias 32:63



ERF8

Pinout of different elements

SAMTEC to Flat converte (by Fabio L.)





ERF8



Pin number SAMTEC	out channel in flat cable
ERF8	
1	N.C.
2	N.C.
3	N.C.
4	N.C.
5	16
6	0
7	17
8	1
9	18
10	2
⋮	⋮
33	30
34	14
35	31
36	15
37	N.C.
38	N.C.
39	N.C.
40	N.C.
41	N.C.
42	N.C.
43	N.C.
44	N.C.
45	48
46	32
47	49
48	33
49	50
50	34
⋮	⋮
71	61
72	45
73	62
74	46
75	63
76	47
77	N.C.
78	N.C.
79	N.C.
80	N.C.

First pin in SAMTEC connectors and cables

ERF8



ERCD



Electrical connection between Bias supply and Tracker

